

SURMODICS VASCULAR DEVICE SOLUTIONS

.035" High Pressure Large Diameter Balloon Dilation Catheters

.035"
HP LD

**HIGH PRESSURE
YOU NEED ... AND A
WHOLE LOT MORE**

LARGE VESSEL DILATATION
WITHOUT COMPROMISE





SURMODICS

Resistant large vessel lesions have met their match!


- Delivering pressures up to 18atm rated burst pressure
- Short balloon shoulders and catheter tip reduce "foot print" in vessel
- Rapid balloon deflation shortens time to post-dilatation imaging

"Vessel-friendly" catheter delivers and treats with ease

- Flexible catheter facilitates tracking through tortuosity
- Conformable balloon lessens vessel straightening during inflation



Product Features



Product Specifications

Catheter Platform	OTW
Recommended Guide Wire	.035"
Sheath Compatibility	7Fr (12mm), 8Fr (14mm & 16mm)
Catheter Working Lengths	75cm
Balloon Diameters	12mm, 14mm, 16mm
Balloon Lengths	40mm
Nominal Pressure	12atm (12mm), 10atm (14mm & 16mm)
Rated Burst Pressure	18atm (12mm), 16atm (14mm), 14atm (16mm)



Creagh Medical and NorMedix are now part of **Surmodics, Inc.**

9924 West 74th Street
Eden Prairie, MN 55344 USA
Toll Free 1-866-787-6639
Phone 952-500-7000
Fax 952-500-7001
www.surmodics.com

© 2016 Surmodics, Inc.
All rights reserved. 35HP0416 01

FOR MORE INFORMATION ON HOW **SURMODICS** CAN ENHANCE YOUR BUSINESS, PLEASE CONTACT US AT **952-500-7000** OR VISIT **WWW.SURMODICS.COM**.